

### Travel to New Dimensions- LSM 700

**ZEISS**

Sensitivity,  
Flexibility and  
Ease of Use

Innovative High-End  
Laser Scanning Microscopes  
from Carl Zeiss

余雅倩  
台灣儀器行



Carl Zeiss Microimaging GmbH, Jena, Germany 2015/1/14 1

### Spatial Resolution of a Light Microscope

**ZEISS**

Objective and tube lens do not image a point as a bright disk with sharply defined edges, but as a slightly blurred spot surrounded by **diffraction rings**

Point Spread Function

$2D_{XY}$


Airy disks

$2D_{XZ}$

Tube lens

Objective

Stage




Carl Zeiss Microimaging GmbH, Jena, Germany 2015/1/14 4

### Basic principle of light microscope

#### Different types of light microscopes

**ZEISS**

Upright Microscope      Inverted Microscope      Stereo Microscope




Carl Zeiss Microimaging GmbH, Jena, Germany 2015/1/14 5

### Basic principle of light microscope

#### Upright Microscope

**ZEISS**

- Short working distance
- Higher magnification
- Higher resolution
- Suitable for slide




Carl Zeiss Microimaging GmbH, Jena, Germany 2015/1/14 6

### Different Beam Path of Image Formation

#### Transmitted-light and Reflected-light in upright microscope

**ZEISS**




Carl Zeiss Microimaging GmbH, Jena, Germany 2015/1/14 7

### Basic principle of light microscope

#### Inverted Microscope

**ZEISS**

- Long working distance
- Cell incubation
- Micromanipulation
- Suitable for petri dish sample



Carl Zeiss Microimaging GmbH, Jena, Germany 2015/1/14 8

**Different Beam Path of Image Formation**  
**Transmitted-light and Reflected-light in inverted microscope**



**Transmitted-light**

**Reflected-light**

Carl Zeiss MicroImaging GmbH, Vanessa

2015/1/14

9

**Basic principle of light microscope**



Carl Zeiss MicroImaging GmbH, Vanessa

2015/1/14

10

**The Cell Observer – New Incubation Concept**  
**Incubation Devices**

**Incubator PM S1**

- Compatible with **petri dishes & multiwell plates**
- Heating inserts: Petri dishes, 6/96 well plates
- Heatable glass surface
- Improved gas mixture delivery
- Fast changes of incubation conditions (in dynamic experiments)
- Compatible with cooling experiments
- Space-saving solution




Carl Zeiss MicroImaging GmbH, Vanessa

2015/1/14

11

**Major tasks of a LSM**  
**Time lapse image**



Carl Zeiss MicroImaging GmbH, Vanessa

2015/1/14


12

**The Point-Spread-Function is a 3-dimensional function**

The axial shape of the PSF is completely different from the lateral one.

The axial extension is larger than the lateral.

→ A microscope has a lateral and an axial resolution.



Carl Zeiss MicroImaging GmbH, Vanessa

2015/1/14

31


**The Resolution of a Microscope is limited**

What does that mean?

The image of a point-like structure is not a point, but a **diffraction pattern** with a finite extension.

This 2-dimensional pattern in the image plane is also called the **Airy-disc**.

In general, the image of a point-like structure is called the **Point Spread Function (PSF)**.



Carl Zeiss MicroImaging GmbH, Vanessa

2015/1/14

32

### The Resolution of a Microscope is limited

**Definition**

The resolution limit is reached, when two point-like objects can not be imaged as two distinct structures anymore.

The distance between the objects is called the resolution limit.

Object

$d=0.4 \mu\text{m}$

Image

Object

$d=0.3 \mu\text{m}$

Image

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 33

### The Resolution of a Microscope is limited

Prof. Ernst Abbe (1840 - 1905)

$d = \frac{\lambda}{2 \cdot NA}$  (1876)

$$d_0 = \frac{1,22\lambda}{NA_{obj} + NA_{cond}}$$

or more simply  $d_0 = \frac{\lambda}{2NA}$

$\lambda$  = wavelength of light, e.g. 550 nm (green)

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 34

### The Resolution of a Microscope depends on Numerical Aperture and Wavelength

$d = \frac{\lambda}{2 \cdot NA}$

$NA = n \cdot \sin \alpha$

**Principle effect**  
High numerical aperture objectives have a large opening angle.

The higher the NA, the better the resolution of the microscope.

Object

$0.3 \mu\text{m}$

Image

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 36

### The Resolution of a Microscope depends on Numerical Aperture and Wavelength

$d = \frac{\lambda}{2 \cdot NA}$

$NA = n \cdot \sin \alpha$

**Principle effect**  
Shorter wavelengths generate smaller Point Spread Functions.

The shorter the wavelength, the better the resolution.

Object

$0.3 \mu\text{m}$

Image

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 38

### Basic principle of light microscope

#### The Numerical Aperture of Objective

The numerical aperture of a microscope objective is a measure of its ability to gather light and resolve fine specimen detail at a fixed object distance.

A. Low Magnification (10X/0.25)

B. High Magnification (40X/0.75)

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 39

### Conventional/Widefield Fluorescence

Background emission from deeper image planes

Structures which are „out-of-focus“ become visible in conventional widefield-fluorescence. Because of the focal depth inherent in all objectives, they are visible as an image blur (haze, image fog).

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 40

### Confocal Principle

Widefield                      Confocal

Carl Zeiss Microscopy GmbH, Dr. Detlef Koch, Training Application and Support Center - TASC APAC Singapore

### General Optical Sectioning Methods

#### General Overview

Carl Zeiss MicroImaging GmbH, Vanessa                      2015/1/14                      42

### Imaging of 3-dimensional objects

#### The fundamental problem

**Conventional Images**

Conventional images of 3-dimensional objects always contain light from structures, which are in focus and light from structures which are not in focus.

Carl Zeiss MicroImaging GmbH, Vanessa                      2015/1/14                      43

### Imaging of 3-dimensional objects

#### The fundamental problem

**Conventional Images**

Conventional images of 3-dimensional objects always contain light from structures, which are in focus and light from structures which are not in focus.

This out-of-focus light blurs the structures from the focal plane and reduces the contrast and resolution.

Carl Zeiss MicroImaging GmbH, Vanessa                      2015/1/14                      44

### The confocal principle

A minute diaphragm, situated in a conjugated focal plane, prevents out of focus light to be detected.

The pinhole diameter directly controls the thickness of the optical section.

Carl Zeiss MicroImaging GmbH, Vanessa                      2015/1/14                      47

### Confocal Laser Scanning Microscopy

#### Optical sectioning: elimination of out-of-focus light

Carl Zeiss MicroImaging GmbH, Vanessa                      2015/1/14                      48

### The Comparison Between the LSM and the Conventional Light Microscope

	Wide Field Microscope	Laser Scanning Microscope
Light Source	Mercury or Xenon Lamp	Laser
Illuminated Field	Wide Field	Spot
Image Acquisition	Parallel, Frame at Once	Sequential, Pixel wise
Signal Separation	Dichroic Beam Splitter, Emission Filter	Beam Splitter Cascade, Emission Filter
Detector	Eye or CCD Camera	Diffraction limited by pinhole → Photomultiplier (PMT)

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 49

### Multiple staining - the crosstalk problem

Simultaneous scan

Simultaneous scan

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 50

### Emission Crosstalk - way around with Sequential image acquisition

Sequential scan

Sequential scan

Vanessa\_Yu Taiwan Instrument Company 2015/1/14 51

### Confocal: Point Scanning

From Spot to Image

- To get a 2 dimensional image from the specimen, the excitation spot has to be moved over the specimen
- The scanning mirrors move the excitation beam in a line wise fashion

XY scanning

Point scanning confocal systems

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 52

### Confocal: Point Scanning Sequential image acquisition

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 53

### Major tasks of a LSM

#### Colocalization in Confocal Microscopy

- Acquisition of **Crosstalk free** images required
- Occurrence of two fluorescent emission signals inside the same detection volume
- Identical size of detection volumes for different color channels required
- Intensities and position of the signals inside the detection volume may vary

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 55

Select all 1 AU or 調整pinhole至相同的光學切片厚度

488 nm

Pinhole 1.00 Airy Units = 0.8 µm section

1 AU max

Carl Zeiss MicroImaging GmbH, Jena 2015/1/14 56

### Major tasks of a LSM

#### Optimal optical sectioning in thick tissue Z stack

This plane represents an optical section

Z-Drive

X/Y/Z Stack

- 3 D information is acquired by moving the excitation focus not only in XY direction but also in Z direction
- The result is a **3 D data stack** consisting of number of XY images representing different optical sections from the specimen

Carl Zeiss MicroImaging GmbH, Jena 2015/1/14 57

### Major tasks of a LSM

#### Optimal optical sectioning in thick tissue Z stack

#### Number of sections

Missing Information

Sample bleached and much data.

Carl Zeiss MicroImaging GmbH, Jena 2015/1/14 58

### Major tasks of a LSM

#### Optimal optical sectioning in thick tissue Z stack

Optimal Number of sections : no missing information at minimal number of sections

Optical thickness depends on:

- wavelength  $\lambda$
- objective lens, N.A.
- refractive index  $n$
- pinhole diameter  $P$

$$d \sim P n \lambda / (N.A.)^2$$

„Nyquist-“ or Sampling- Theorem  
slices overlap by the 50% of their thickness

LSM software: One click for best resolution

Carl Zeiss MicroImaging GmbH, Jena 2015/1/14 59

### Major tasks of a LSM

#### Optimal optical sectioning in thick tissue

0 µm 2 µm 4 µm 6 µm 8 µm  
10 µm 12 µm 14 µm 16 µm 18 µm  
20 µm 22 µm 24 µm 26 µm 28 µm

A series of of confocal images from different optical planes contains the image information from the whole specimen

Carl Zeiss MicroImaging GmbH, Jena 2015/1/14 60

### Major tasks of a LSM

#### Optimal optical sectioning in thick tissue

- An overlay (maximum projection) of these single images results in an image with an **enhanced depth of focus**
- This image contains all information from the specimen

↓

**Every detail is in focus !**



Carl Zeiss MicroImaging GmbH, Jena 2015/1/14 61



### The Power of Sensitivity

Brighter Images ≠ Increased Sensitivity !

Signal to noise ratio is critical !

$$SNR = \frac{S}{N}$$





Carl Zeiss Microimaging GmbH, Jena  
2015/1/14 62

### The Power of Sensitivity

Brighter Images ≠ Increased Sensitivity !

Signal to noise ratio is critical !

An increase in brightness of the image does not provide better information

Sources of noise:

- Shot noise

Sources of noise:

- Dark noise
- Electronic noise
- Laser reflection



Carl Zeiss Microimaging GmbH, Jena  
2015/1/14 63

### The Power of Sensitivity

Brighter Images ≠ Increased Sensitivity !

Signal to noise ratio is critical !

A decrease in noise gives cleaner images



Carl Zeiss Microimaging GmbH, Jena  
2015/1/14 64

### The Power of Sensitivity

Brighter Images ≠ Increased Sensitivity !

Signal to noise ratio is critical !

Low-noise images can be displayed with optimal brightness


Noise reduction in LSM 700:

- Reduced dark noise
- Reduced electronic noise
- Reduced laser reflection

Carl Zeiss Microimaging GmbH, Jena  
2015/1/14 65

### LSM 700 in Standard Applications

- Uncompromised Image quality



Improved signal recording:  
Crisp details, clear image data

Improved S/N ratio:  
Black background

Carl Zeiss Microimaging GmbH, Jena  
2015/1/14 66

### LSM 700 Laser line


Laser line	Fluorochrome
405 nm	DAPI, Hoechst, Alexa 405, BFP,
488nm	Alexa 488, Fluo-4, FITC, eGFP
555 nm	Rhodamine, Alexa 546, 555, 568, Cy3, TRITC, DsRed, Texas Red, MitoTracker Red
639 nm	Alexa 633, Cy5..

Detectors:  
2 reflection PMT detectors for fluorescence images  
1 transmitted PMT detector for Bright Field(PH/DIC) images

Carl Zeiss Microimaging GmbH, Jena  
2015/1/14 68


### LSM 700


#### VSD – Variable Secondary Dichroic



- VSD is a variable short pass beam splitter for splitting signals between detectors
- Positioning of VSD allows precise tuning of wavelength at which signals are split (splitting possible between 420 and 630 nm, min. step: 1 nm)
- Enables highly light-efficient detection strategies and spectral imaging (lambda stack acquisition)
- Patented Zeiss innovation

**Principle of the VSD**







Carl Zeiss Microimaging GmbH, Jena 2015/1/14 69

### LSM 700

#### VSD – the new flexible way





**Example: Dual-color Detection of GFP and MitoTracker Orange**

**Approach: "New Flexible Way"**


Flexible dual-color detection enabled by the new variable secondary dichroic (VSD) of the LSM 700.


**Improvement: Enhanced light efficiency because no portion of the signal is excluded from the detection process.**

Carl Zeiss Microimaging GmbH, Jena 2015/1/14 74

### LSM 700

#### VSD – the new flexible way





**Example: Dual-color Detection of GFP and MitoTracker Orange**

**Approach: "New Flexible Way"**


Flexible dual-color detection enabled by the new variable secondary dichroic (VSD) of the LSM 700.


**Improvement: Enhanced light efficiency because no portion of the signal is excluded from the detection process.**

Carl Zeiss Microimaging GmbH, Jena 2015/1/14 75

### LSM 700

#### VSD – the new flexible way





**Example: Dual-color Detection of GFP and MitoTracker Orange**

**Approach: "New Flexible Way"**



Flexible dual-color detection enabled by the new variable secondary dichroic (VSD) of the LSM 700.

**Also possible: use of emission filters (optional) for additional specificity.**

Carl Zeiss Microimaging GmbH, Jena 2015/1/14 76

### ZEN 2011 - Efficient Navigation


Powerful software for powerful LSM systems






Carl Zeiss Microimaging GmbH, Jena 2015/1/14 87

### ZEN 2011

#### Load configuration







Carl Zeiss Microimaging GmbH, Jena 2015/1/14 88



### ZEN 2011 Load configuration

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 89

### Major tasks of a LSM Laser and scanning mirror control

- Easy sample manipulation
- Flexible scanning strategies (1D to multiD)  
Scan Mode --- 1D, 2D, and free 2D Image

Absolut linear scanner movement:  
The same dwell time for every pixel in the images  
(essential for any quantitative measurements)

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 90

### Major tasks of a LSM Laser and scanning mirror control

**Real Regions of Interest (rROI)**  
Irregular shaped areas  
Up to 99 areas simultaneously  
Sample irradiation only during data Acquisition (beam blanking)  
No photobleaching in surrounding areas

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 91

### Major tasks of a LSM Laser and scanning mirror control

**Two independent scanning mirrors**

Free scan field rotation (0-360°)  
Free online zooming (0.6-40x (zoom=66.7x))  
Any geometry: 1x4... 6144\*6144  
Faster rectangular acquisition (e.g. video rate)

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 92

### Tile scanning with motorized scanning stage 大面積高倍數掃描

40X objective, 10X9

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 93

### Major tasks of a LSM Laser and scanning mirror control Photomanipulation for studying cellular dynamics

**Photomanipulation**

- Photobleaching
- Photoactivation
- Photoconversion
- Uncaging
- Laser Ablation

Carl Zeiss Microimaging GmbH, Vanessa 2015/1/14 94

### Photoconversion - Kaede

**An optical marker based on the UV-induced green-to-red photoconversion of a fluorescent protein**  
 Ryoko Ando<sup>1</sup>, Hiroshi Hamai<sup>2</sup>, Miki Yamamoto-Rino<sup>2</sup>, Hirosaki Mizuno<sup>2</sup>, and Atsushi Miyawaki<sup>1\*</sup>  
 PNAS | October 1, 2002 | vol. 99 | no. 20 | 12651-12656

Kaede = maple leaf (jap.)  
 New fluorescent protein from the stony coral *Trachyphylla geoffrey*  
 Includes Tripeptide which acts as green chromophore that can be converted to red  
 Conversion highly sensitive to irradiation with UV light  
 Example: allows to delineate a single neuron in a dense culture

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 95

### Photoconversion - Kaede

**Photoconversion from green to red**

**Sample**  
 Cell culture stable expressing cytoplasmatic KAEDE

Convertible Fluorescent Protein (from green to red)

Recorded time: 45 sec

Conversion: 405 Laser  
 Left region: 50%  
 Right region: 100%

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 96

### 1. 開關機步驟

1. 總電源(延長線上)

2. 螢光燈源

3. 電腦電源

開機1→4順序開啟  
 關機4→1順序關閉

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 97

### 1. 進入軟體ZEN

### 2. 開啟軟體與硬體連結

進階模式

啟動機器掃圖

\*啟動軟體全功能，但不與硬體做連結，單純分析資料、沒有要操作機器請選此項。

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 98

### 3. 控制顯微鏡找到樣品焦距

1. Online，自目鏡觀察，找到樣品焦距。  
 2. 利用configuration切換各贏光濾片組與穿透光設定。  
 3. 確認樣品位置及焦距後切至Offline，即可以進入LSM影像擷取模式。

Online: 分光至目鏡  
 Offline: 分光至LSM，此時目鏡無法做觀測

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 99

### 4. Apply Configuration Setting

套用適合的Configuration (選擇染劑名字)

Carl Zeiss Microimaging GmbH, Jena, 2015/1/14 100

### 5. Acquisition setting

Select all  
1 AU或調整pinhole至相同的光學切片厚度

488 nm  
Pinhole  
1.00 Airy Units = 0.8 μm section  
1 AU max

2015/1/14 101

### 5.1 :設定適當的Pinhole大小

488 nm  
Pinhole  
1.00 Airy Units = 0.8 μm section  
1 AU max

Pinhole設定:  
一般選擇1AU  
亦可以設定成所有channel 為相同的um厚度, 例如在高倍物鏡下可以統一選擇1um。

2015/1/14 102

### 6. 預覽掃圖

1 Snap 試拍一張  
2 Continuous  
按**Continuous**後進入預覽模式, 同時使用focus drive精調螢幕中影像焦距  
3 STOP  
找到焦距後停止掃描

2015/1/14 103

### 6.1 掃描參數值設定

1 & 2. 於frame size 512\*512畫素、speed 8下進行影像快速掃描, 以方便即時預覽更改參數後的結果

- 3. Continuous: 持續掃描, 要按停止才會停止掃描
- 4. Snap: 拍一張影像

2015/1/14 104

### 6.2 掃描參數值設定

- 依序設定每個track的掃描參數
- 調整PMT gain值, 數值越大訊號被放大得越多, 影像越亮, 可使用滑鼠中鍵滾輪滑動調整。  
(建議數值600-750)
- 調整laser強度, 數值越大, 影像越亮
- 調整好後進行下一個track的設定, 重複1~3步驟直到每個track都設定完畢。



2015/1/14 105

### 6.3 掃描參數值設定

拍圖調整掃描條件時建議選取range indicator套色方式表現色彩, 將有助於將顏色之intensity調到最佳分布。紅色表示飽和(調整detector gain和laser量), background 深藍表示全黑(調整digital offset)  
建議調整到全畫面當中訊號少部分飽和, background 部分為藍色。

2015/1/14 106

### 6.4 掃描參數值設定 ZOOM、掃描區域選擇



- Corp功能包含zoom in/ zoom out
- Zoom 勿過度使用,一般不會超過3 否則將造成bleach樣品的效果

Carl Zeiss Microimaging GmbH, Venzona

2015/1/14

107

### 7. 正式掃圖:獲得高解析度影像



1. 將所有已經調整好參數的track都打勾
2. 選擇需要的畫素,一般需要發表須要1024\*1024.
3. 調整掃圖至慢速度,高品質影像建議scan speed為5~7
4. 一般均設定1,若影像品質不佳可採用平均數次可以使影像品質提升,降低雜訊
5. SNAP拍一張,獲得漂亮的data!!




Carl Zeiss Microimaging GmbH, Venzona

2015/1/14

108

### 8. Z stack



1. 進入Z stack
2. 進入 Mark First/Last
3. Continuous scanning
4. 用 Focus drive找到觀察樣品厚度之最高/低點→mark first
5. 反方向轉動粗細調節輪找到欲觀察範圍的最高/低點→mark last
6. 滑鼠按下 Optimal interval建議值
7. 設定完畢後stop,避免樣品被bleach
8. 回到掃圖設定成1024\*1024, speed 7-5
9. Start experiment 開始執行Z section拍照



Carl Zeiss Microimaging GmbH, Venzona

2015/1/14

109

### 9. 存檔—\*.lsm完整檔案,可以reuse!

File → Export




Format : 下拉選擇欲存檔之類型 (tif、jpg...), 建議以tif儲儲存,減少壓縮損失,亦可儲存serial section 成動畫影片檔(.avi)

Carl Zeiss Microimaging GmbH, Venzona

2015/1/14

110

### 10. 輸出成圖片檔或者影片檔



- Raw data : 不含尺規,選擇要存的顏色,是否為灰階等等。
- Content of image window : 存下室窗內的影像畫面,包含尺規。
- Full resolution : 包含尺規依照拍照時的畫素存檔 (建議使用!)

- 建議使用Full resolution 或者Contents of image window
- single plane : 單張,目前所顯示的單層/單張影像。
- series : 一系列圖,適用於Z stack, time series和movie檔。


Carl Zeiss Microimaging GmbH, Venzona

2015/1/14

111

### 11. 把多張Z section疊成一張

製造全景深影像: maximum intensity projection



Carl Zeiss Microimaging GmbH, Venzona

2015/1/14

112