

16S Analysis on BaseSpace




Jora Lin 林孟儒
Bioinformatic Specialist
均泰生物科技有限公司

techsupport@gtbiotech.com.tw

Demultiplexing

Multiplexed samples are sequenced together
Unique index sequences allow sequencing data reads to be sorted by samples



What is The BaseSpace?

BaseSpace is a website in the 'Cloud'



- **BaseSpace** is hosted in **Amazon** Web Services
- Anyone can use **BaseSpace** – just go to basespace.illumina.com
 - Available for **illumina** customers
- **BaseSpace** currently supports
1TB free storage for each account

The Benefit for BaseSpace

▶ Online Monitor

✓ Check your sequencing run at anytime with Internet

▶ Backup

✓ Automatic Online Backup

▶ Analysis App



✓ Very easy to use for anybody


▶ Data Sharing

✓ No HDD or USB storage, just share data by account

▶ Demo Data


✓ The Illumina data updated






Remote Instrument Run Monitoring

Requires Web Connection





Sharing and Collaboration

Transfer Ownership
 ■ NewProject

New Owner

jsmith@email.com

Optional Message

Hello Jane,

I am transferring ownership of this project and associated analysis to you. After the transfer I will have no access to the project

ⓘ Apps in your account that currently use this data may not function properly after the data is transferred out of your account. Please review the [BaseSpace User Terms of Use](#) for more information.

Cancel Continue

Share This Project
 ■ NewProject

Invite A Collaborator

jsmith@email.com

Hi Jane,
 I would like to share the analysis from this project with you

Add Collaborator

Collaborator Settings

jsmith@email.com
 pending

Read Only

Cancel Save Settings

Share By Link
 ■ NewProject

Copy the below URL and send it to anyone you want to share with.

https://basespace.illumina.com/s/4TpIG0WdFqqf

Close Deactivate

Easily share or transfer ownership of runs, samples, or projects with anyone in the world*

*Sequence Hub account with available storage required to accept transfer

3



Push-Button Analysis

Three types of pre-built Apps to analyze your data

The screenshot shows the BaseSpace interface with the 'APPS' tab selected. Under the 'Applications' heading, three categories are displayed:


- Core Apps:** Fully Supported, Illumina Developed Apps. Example: Cufflinks Assembly & DE (Illumina, Inc.).
- Labs Apps:** Apps created by Illumina R&D Groups. Example: Integrative Genomics Viewer (BaseSpace Labs).
- 3rd Party and Custom Apps:** See App page for support. Example: Genomatix Pathway System (GePS) (Genomatix).

At the bottom of the interface, a green banner states: **Over 80 published Apps supporting all of Illumina library prep kits**




Push-Button Analysis

Common Workflow for All Illumina Core Apps



Analysis Info	
Name	Whole Genome Sequencing 03/21/2017 5:19:53
Application	Whole Genome Sequencing Version: 4.0.0
Date Started	Tuesday, March 21 2017 5:24:50 PM
Date Completed	N/A
Duration	6 minutes 6 seconds
Compute Charge	0.00 iCredits
Session Type	Multi-Node
Status	Running (1 Node Pending Execution)



- Sequencing Summary Report
- Alignment Map
- *.bai (Index for Alignment Map)
- Variant Call Files
- Log and Error Files
- Summary Reports

BaseSpace performs alignment, calls and annotates variants, summarizes results

BaseSpace Public Data

Public Data [?](#)

[iSeq 100: Nextera DNA Flex \(E.coli, B.cereus, Rhodobacter\)](#)

Metagenomics

[iSeq 100: Cancer Hotspot v2 Panel \(Replicates of NA12878\)](#)

targeted-sequencing variant-analysis ampliseq

[HiSeq 2500](#)

resequencing

MiSeq: Nextera DNA Flex (replicates of E.coli, B.cereus, and R.sphaeroides)

3 replicates each of B. cereus, E. coli, and R. sphaeroides at 1ng and 100ng DNA input were prepared using the Nextera DNA Flex library prep kit and sequenced at 2x151 on a MiSeq instrument.

resequencing denovo-assembly

Run	MiSeq_Nextera_DNA_Flex_E-coli_B-cereus_R-sphaeroides (12 GB)	Import
Project	MiSeq: Nextera DNA Flex (replicates of E.coli, B.cereus, and R.sphaeroides) (10.83 GB)	Import

[NovaSeq S4: IDT xGen Exome Research Panel \(96 replicates of NA12878\) – UDI](#)

resequencing


[NovaSeq S1 Xp: TruSeq Nano 350 \(Replicates of NA12878\)](#)

resequencing

Categories

Exome (12)	Resequencing (49)
Tumor Normal (1)	Variant Analysis (15)
Differential Expression (5)	Quality (0)

Register a BaseSpace account



BaseSpace **SEQUENCE HUB** illumina

BaseSpace Sequence Hub is a cloud-based genomics analysis and storage platform that directly integrates with all Illumina sequencers.


[Learn More](#)

[Sign up](#) [Log in](#)

- Set up runs for Illumina sequencers and NeoPrep
- Monitor sequencing runs from the web or with the iOS app
- Stream data to the cloud directly from sequencers
- Analyze genomics data with pre-configured pipelines
- Share data instantaneously with anyone in the world

[Developers](#) [Blog](#)

basespace.illumina.com



illumina

Register

Already have an account? [Sign in](#)

Email address


First name

Last name

Password [Guidelines](#)

Confirm password

Country
 United States

I'm not a robot  [Privacy - Terms](#)

By creating an account, you agree to Illumina's [Terms of Service](#)

[Create Account](#) [Cancel](#)

BaseSpace Dashboard

BaseSpace

SEQUENCE HUB

DASHBOARD

MY DATA

PREP

APPS

PUBLIC DATA

Dashboard Personal

Your FREE TRIAL has expired

Basic Tier

Usage

iCredits subscription required

Ready for more? Upgrade to a Professional subscription.

[Learn More](#)

Storage [Details](#)

Personal Storage in use **14.31 GB**

Developers [Details](#)

Start building the next generation of BaseSpace sequencing apps.

Newsfeed

BaseSpace Cohort Analyzer Update
BaseSpace Cohort Analyzer enables users to apply complex genomic data in novel...
Feb 24, 2017

Updated Command Line Tools - Wait for App Dependencies
We are pleased to announce a minor release of BaseSpaceCLI (0.8.10) with...
Dec 6, 2016

BaseSpace Informatics Suite Summit 2016
You're invited to an exclusive informatics event Advancing Precision...
Sep 10, 2016

VCAT v3.0 and FASTQ Toolkit v2.2 - Health Improvements and Gold

Notifications

- Workgroups** NEW
Mar 30, 2017
Workgroup users now see their last-viewed workgroup when signing in to BaseSpace Sequence Hub.
- Share accepted**
Jan 25, 2017
You accepted Dorothea Aglus's invitation to the Project MiniSeq: TruSight...
- Share accepted**
Jan 23, 2017
You accepted Illumina Inc's invitation to the Project MiSeq: TruSight RN...
- Share accepted**
Jan 17, 2017
You accepted Eric Allen's invitation to the Project MiSeq v3: 16S Meta...
- Share accepted**
Jan 11, 2017
You accepted BaseSpace User2's invitation to the Project Illumina FastTrack...

Latest Runs

0 None at this time

Your latest sequencing runs and their statuses will appear here.

[All Runs >>](#)

Latest Analyses

- Complete**
Jan 25, 2017
20161217MG02_S2
- Complete**
Jan 23, 2017
RNA-Seq Alignment 01/23/2017 12:01:52
- Complete**
Jan 23, 2017
RNA-Seq Alignment 01/23/2017 12:06:22
- Aborted**
Jan 20, 2017
ChIPSeq 01/11/2017 11:06:16
- Complete**
Jan 18, 2017
16S Metagenomics 01/17/2017 6:34:25

[All Analyses >>](#)

Report

Analysis Summary and result of each sample

Projects

Biosamples

Runs

Analyses

MISeq v3: TruSight Tumor 15 - Sample HDC749

Analysis: TruSight Tumor 15 04/16/2018 1:17:34



Launch App

File

Status

Share

Move to Trash

SUMMARY **REPORTS** INPUTS FILES

Summary

Pair 6: HDC749_Rep6 Report

Pair 1: HDC749_Rep1 Report

Pair 3: HDC749_Rep3 Report

Pair 2: HDC749_Rep2 Report

Analysis Summary Report

Pair 4: HDC749_Rep4 Report

Pair 8: HDC749_Rep8 Report

Pair 7: HDC749_Rep7 Report

Pair 5: HDC749_Rep5 Report

RESULTS FOR PAIR 6: HDC749_REP6

PDF Summary Report

HDC749_Rep6_FilteredVariantSumm... 1 / 2

Sample ID	HDC749_Rep6	Analysis Date	2018-04-16 21:12
Run Name	N/A	Report Definition	TST_15-ReportDefinition-v2.0.0

The Filtered Variant Summary is designed to list a specific subset of variants in the genes targeted by TruSight Tumor 15. A list of these pre-defined variants can be found in the Report Definitions file. Refer to the sample's genome.vcf and amplification.vcf files for an unfiltered list of all variants found in the regions targeted by TruSight Tumor 15. The Report Definitions file and sample vcf files are found in the Analysis Output folder.

Analysis Results

Detected Variants

GENE	VARIANT OF INTEREST	COORDINATE	NUCLEOTIDE CHANGE	CONSEQUENCE	FREQUENCY	COVERAGE
<i>BRAF</i>	p.V600E	chr7:140453136	c.1799T>A	missense variant	0.122	17060
<i>EGFR</i>	p.G719S	chr7:55241707	c.2155G>A	missense variant	0.266	14938
	p.L858R	chr7:55259515	c.2573T>G	missense variant	0.039	845
<i>KRAS</i>	p.G12D	chr12:25398284	c.35G>A	missense variant	0.066	3861
	p.G13D	chr12:25398281	c.38G>A	missense variant	0.167	3863
<i>NRAS</i>	p.Q61K	chr1:115256530	c.181C>A	missense variant	0.096	6797
<i>PIK3CA</i>	p.E545K	chr3:178936091	c.1633G>A	missense variant	0.075	5337
	p.H1047R	chr3:178952085	c.3140A>G	missense variant	0.208	10267

Output Files

Results organized into sample folders

Projects

Biosamples

Runs

Analyses

MISeq v3: TruSight Tumor 15 - Sample HDC749

Analysis: TruSight Tumor 15 04/16/2018 1:17:34

SUMMARY

REPORTS

INPUTS

FILES

Launch App




File

Status

Share

Move to Trash

Files

NAME	ITEM(S) ANALYZED	DATE CREATED	TYPE	SIZE
 Pair 8: HDC749_Rep8	2 items	4 hours ago	App Result	161.06 MB
 Pair 7: HDC749_Rep7	2 items	4 hours ago	App Result	202.55 MB
 Pair 6: HDC749_Rep6	2 items	4 hours ago	App Result	190.33 MB





























bam, VCF file...

BaseSpace App

BaseSpace **SEQUENCE HUB** DASHBOARD MY DATA ▾ PREP **APPS** PUBLIC DATA 🔍 ? 🗑️ 👤


Search Apps 🔍 All Categories ▾

All Sequence Hub Applications


 16S Metagenomics Illumina, Inc.	 Amplicon DS Illumina, Inc.	 Bacterial Analysis ... GoSeqIt ApS	 BWA Aligner BaseSpace Labs
 BWA Enrichment Illumina, Inc.	 Cancer Variant Cal... Samsung SDS	 ChIPSeq BaseSpace Labs	 CrossMap BaseSpace Labs
 Cufflinks Assembl... Illumina, Inc.	 DeepChek@-HBV ABL S.A.	 DeepChek@-HCV ABL S.A.	 DeepChek@-HIV ABL S.A.
 DESeq2 BaseSpace Labs	 DNA Amplicon Illumina, Inc.	 DRAGEN Germline ... Edico Genome Inc.	 E. coli Serotyping GoSeqIt
 EDGC Annotator EONE-DIAGNOMICS G...	 Elastic Genome Br... Strand Life Sciences	 Enrichment Illumina, Inc.	 FASTA Upload BioBam Bioinformatic...
 FASTQ Toolkit BaseSpace Labs	 FastQC BaseSpace Labs	 GeneTalk Variant A... GeneTalk GmbH, Berlin	 GENIUS Metageno... CosmosID
 GePS Genomatix Pathwa... Genomatix	 Genome Profiler Biomatters Ltd	 Hap.py Benchmark... BaseSpace Labs	 IgReC Center for Algorithmic...

contact us

16s Metagenomics Analysis



16S rRNA sequencing with Illumina protocols & software tools



BaseSpace

16S Metagenomics App



- ▶ GreenGenes database to perform taxonomic classification
 - Illumina-curated version
- ▶ Bayesian classification method to assign taxonomies
 - Short sub-sequences (32 bases) are extracted from each read and compared to the database by the classifier
 - Classification down to genus/species-level

Demultiplexing RAW Files, FASTQ File

Instrument Name

Run number

Flowcell ID

Lane

Tile

X_pos

Y_pos

Read

Is Filtered

Control Number

Index Sequence

@HWI-BRUNOP20X:994:B809UWABXX:1:1101:13501:2240 1:N:0:CTTGTA

TGAAACCAGTGTTCCTTAATTGGCATTTCACACACACACACAGAATTTAAAAAAAAAATCAAAGGAAATCATTCTAAATGTACTATGATAGCATGTTAAA

+

=55>7;?::BDADDD@EE88DCD?DFFEFFECBE6666BB=B;<;<-34:;<CB51>=BBEE>EE?3D@??CB->=:AA8DDDDDBBE9; ,=?:/89<E

$$Q\text{-score} = -10 \log(\text{error rate})$$

Phred Score

+ 33

ASCII value


Q-score=Phred Score

Character	ASCII Value	Phred Score	Error probability
5	53	20	0.01
?	63	30	0.001
!	73	40	0.0001

Examples of ASCII in Sanger FastQ

Read Classification of 16s Metagenomics app

- The original RDP classifier algorithm used 8-base words due to implementation constraints. ClassifyReads uses more efficient data structures and is able to use 32-base words – giving each word more specificity for each species.



Reference sequence



Appendix – Filters applied to GG by Illumina

1. Filter all entries where the 16S sequence length was below 1250 bp
2. Filter all entries that had more than 50 wobble bases
(i.e. M, R, W, S, Y, K, V, H, D, B, N)
3. Filter all entries that were only partially classified
(no classification for genus or species)



16s Metagenomics app



16S Metagenomics v1.0.1
Illumina, Inc.

Analysis Name: i



Save Results To: Select Project(s): i

20171013_16s analysis on BaseSpace	x
------------------------------------	---

Sample(s): Select Sample(s): i

41	x
44	x
47	x
42	x
48	x

AGGREGATE RESULTS



AGGREGATE RESULTS

- PDF Aggregate Summary Report
- Kingdom-Level Aggregate Results
- Phylum-Level Aggregate Results
- Class-Level Aggregate Results
- Order-Level Aggregate Results
- Family-Level Aggregate Results
- Genus-Level Aggregate Results
- Species-Level Aggregate Results

SAMPLE INFORMATION

Sample Number	Sample ID	Number Reads PF	% Reads PF Classified to Genus
1	41	41,862	74.16%
2	44	50,075	57.13%
3	47	110,434	67.41%
4	42	100,803	64.93%
5	48	67,342	65.76%



- ▶ Shannon species diversity measures the entropy of Species-level classifications in the sample.

Sample Number	Sample ID	Shannon Species Diversity	Number of Species Identified
1	Alum_Rock_Cave	1.734	1,440
2	Human_Saliva	3.074	1,111
3	Wastewater_Sludge	2.799	1,789

http://en.wikipedia.org/wiki/Shannon-Wiener_index



RESULTS FOR SAMPLE

CLASSIFICATION STATISTICS

Taxonomic Level	Reads PF Classified to Taxonomic Level	% Reads PF Classified to Taxonomic Level
Kingdom	109,372	99.04%
Phylum	102,850	93.13%
Class	98,852	89.51%
Order	96,184	87.10%
Family	78,502	71.08%
Genus	74,448	67.41%
Species	51,579	46.71%

TOP 20 CLASSIFICATION RESULTS BY TAXONOMIC LEVEL

This column chart shows the relative abundance of the top 20 classification results within each taxonomic level. Mouse over any category to see its description and abundance.



Output Files

- [Sample.report.pdf](#)
 - Sample Information and Classification Information
 - [Sample.summary.csv](#)
 - summarized counts of how reads were classified for each level in the sample
 - [Sample.txt.gz](#)
 - The classification assigned to each read
 - [Class-level-Aggregate_count.csv](#)
 - Per-level aggregate counts for all samples, each row represents a unique classification that occurred in one or more samples.
 - [MetagenomicsAggregateReport.pdf](#)
 - Sample Classified Information and Species Diversity Results
-

BaseSpace

QIIME App



QIIME Preprocessing
QIIME Development Team

- **OTU picking, taxonomic assignment**
- **sl-out** folder: contain the summary of read quality file – `split_library_log.txt`
- **closed-ref** folder: contain the OTU table – `otu_table.biom` and tree file – `97_otus.tree`
- If your reads will not be well represented by the Greengenes database, that might be underrepresented




QIIME Visualizations
QIIME Development Team

- **diversity analyses and visualizations**
- Beta diversity plots
- Alpha rarefaction plots
- Taxonomic summaries

Read Classification of QIIME app

uclust_ref (Edgar, RC 2010), as uclust, but takes a reference database to use as seeds.
Taxonomy database: Greengenes v13.8


Searching



- ✓ searching a database stored in memory
- ✓ A query sequence matches a database sequence if the identity is high enough
- ✓ The minimum identity is set by the --id option, e.g. --id 0.97

ATCGGAC 5 matches, 1 mismatch (default allow 1 mismatch)
A-CAGAC **Identity= matches(5)/shorter length(6)=83.3%**

Clustering



- ✓ Query sequences are processed in input order. If a match is found to a database sequence, then the query is assigned to that cluster

Output data form QIIME Visualizations



名稱	修改日期
arare_max10000	2017/5/9 下午 03...
bdiv_even10000	2017/5/9 下午 03...
taxa_plots	2017/5/9 下午 03...
biom_table_summary.txt	2017/5/9 上午 11...
index.html	
log_20150806201556.txt	2017/5/9 上午 11...
table_even10000.biom.gz	2017/5/9 上午 11...
table_mc10000.biom.gz	2017/5/9 上午 11...


Run summary data	
Master run log	log_20150806201556.txt
BIOM table statistics	biom_table_summary.txt
Filtered BIOM table (minimum sequence count: 10000)	table_mc10000.biom.gz
rarefied BIOM table (sampling depth: 10000)	table_even10000.biom.gz
Taxonomic summary results	
Taxa summary bar plots	bar_charts.html
Taxa summary area plots	area_charts.html
Alpha diversity results	
Alpha rarefaction plots	rarefaction_plots.html
Beta diversity results (even sampling: 10000)	
PCoA plot (weighted_unifrac)	index.html
Distance matrix (weighted_unifrac)	weighted_unifrac_dm.txt
Principal coordinate matrix (weighted_unifrac)	weighted_unifrac_pc.txt
PCoA plot (unweighted_unifrac)	index.html
Distance matrix (unweighted_unifrac)	unweighted_unifrac_dm.txt
Principal coordinate matrix (unweighted_unifrac)	unweighted_unifrac_pc.txt

官方使用BaseSpace QIIME分析完的結果可於此處下載:


<https://basespace.illumina.com/s/U8W5Qlq1SyfT>

Output data form QIIME Visualizations

Taxa summary bar plots for each level



Taxa summary area plots for each level




Legend		Total	2	4	5	6	7	10	11	13	14	15	18	19	22	23	24	25	26	29	31	33	36	37	41	42	44	47	48	
Taxonomy		%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	%	
■	k__Archaea;p__Crenarchaeota	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
■	k__Archaea;p__Euryarchaeota	0.1%	0.0%	0.1%	0.0%	0.0%	0.1%	0.0%	0.0%	0.2%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.4%	0.0%	0.1%	0.0%	0.1%	0.1%	0.3%	0.2%	0.0%	0.0%	0.0%	0.0%	0.1%
■	k__Archaea;p__[Parvarchaeota]	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
■	k__Bacteria;p__	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
■	k__Bacteria;p__AC1	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
■	k__Bacteria;p__Acidobacteria	0.3%	0.2%	0.4%	0.1%	0.1%	0.1%	0.2%	0.3%	0.6%	0.3%	0.3%	0.1%	0.1%	0.1%	0.1%	0.4%	0.9%	0.2%	0.7%	0.2%	0.6%	0.4%	0.7%	0.4%	0.2%	0.1%	0.1%	0.2%	
■	k__Bacteria;p__Actinobacteria	8.5%	4.8%	18.6%	3.8%	5.1%	19.7%	12.8%	17.8%	17.2%	2.9%	3.3%	7.4%	8.5%	7.2%	7.5%	8.1%	12.8%	6.0%	13.7%	5.1%	10.1%	4.7%	5.8%	7.7%	3.4%	2.5%	6.8%	6.8%	

Output data form QIIME Visualizations

Alpha diversity results


Alpha rarefaction plots



Beta diversity results

PCoA plot (weighted_unifrac)






Control panel ↓







Resource

- QIIME Preprocessing: <https://basespace.illumina.com/s/Tk0L7PDydgQs>
 - QIIME Visualizations: <https://basespace.illumina.com/s/IYcBNizGBcW8>
 - Both: <https://basespace.illumina.com/s/U8W5Qlq1SyfT>
 - QIIME Forum: <http://forum.qiime.org/>
-

DNA base sequencing

Library Prep Name	App Name
<p data-bbox="214 444 585 486">TruSeq Nano DNA</p> <p data-bbox="214 508 672 551">TruSeq DNA PCR-Free</p>	<p data-bbox="1010 372 1136 494"> Whole Genome Sequencing</p> <p data-bbox="1010 508 1136 629"> Tumor Normal</p>
<p data-bbox="214 758 716 801">TruSeq DNA Methylation</p> <p data-bbox="214 822 624 865">TruSeq Methyl EPIC</p>	<p data-bbox="1010 751 1136 872"> MethyIseq</p> <p data-bbox="1363 751 1489 872"> Methylkit</p>
<p data-bbox="214 1093 469 1136">TruSeq ChIP</p>	<p data-bbox="1010 1051 1136 1172"> ChIPSeq</p>

RNA base sequencing

Library Prep Name	App Name
<p>TruSeq RNA</p> <p>TruSeq Stranded mRNA</p> <p>TruSeq Stranded Total RNA</p> <p>TruSeq RNA Access</p>	<p> RNA-seq Alignment</p> <p> Cufflinks Assembly & DE</p>
<p>TruSeq Targeted RNA Expression</p>	<p> TruSeq Targeted RNA</p>
<p>TruSeq Small RNA</p>	<p> Small RNA</p>

Data preprocessing



FASTQ Toolkit



FastQC

Metagenomics



16s Metagenomics



Prokka Genome Annotation



Kraken Metagenomics

De novo assembly



Velvet de novo Assembly

Visualization



Integrative Genomics Viewer

Public database



SRA Import



SRA Submission

Partek Introduction

Jora Lin 林孟儒

Bioinformatic Specialist


均泰生物科技有限公司

techsupport@gtbiotech.com.tw




Who is Partek?


- Founded in 1993
- Building tools for statistics & visualization
- Focused on genomics
- Thousands of customers worldwide
- Worldwide, world-class customer support



Comprehensive Solution for NGS Analysis



- Flow Server
 - Core engine
- Alignment toolkits
 - Different choice of aligners
 - QA/QC on pre- and post-alignments
- Applications (NGS data)
 - RNA-seq, DNA-seq, smallRNA-seq, ChIP-seq, MeDIP-seq
- Customizable protocols
 - Collection of tasks
 - Protocol sharing
- Enable collaboration / project sharing
- Project/Data management
- Accessed from anywhere with internet connection



Partek Flow Queue Projects

Home > Variant Test


Analyses Data Log Project settings Attachments

	Sample name	Files	
		fastq	+/-
1	UHRRmRNA_S5	UHRRmRNA_S5_L001_R1_001	


Hide data files Download

Import data Assign sample attributes from a file

Add a system-wide attribute column Manage attributes




- Built on statistical and visualization package
- Easy to use workflow by providing step-by-step guidance on how to analyze your data
- Support for integrative analysis between different datasets
- Provide integrative analysis with publicly available database.



The screenshot displays the software's workflow selection interface. On the left, a 'Workflows' dropdown menu is open, listing various analysis types: Choose..., Gene Expression, MicroRNA Expression, Exon, Copy Number, Allele Specific Copy Number, LOH, Association, Trio, Tiling, ChIP-Seq, RNA-Seq, miRNA-Seq, DNA-Seq, Methylation, and High Throughput Screening. On the right, a detailed panel for the 'RNA-Seq' workflow is shown. This panel includes sections for 'Import' (with options like 'Import and manage samples', 'Add sample attributes', and 'Choose sample ID column'), 'QA/QC', 'Analyze Known Genes' (with options like 'mRNA quantification', 'Differential expression analysis', 'Alternative splicing analysis', and 'Create gene list'), 'Allele-Specific Analysis' (with 'Detect Single Nucleotide Variations'), 'Visualization' (with 'Plot chromosome view' and 'Cluster based on significant genes'), and 'Biological Interpretation' (with 'Gene set analysis' and 'Pathway analysis').

Partek® Pathway™?

- Find Enriched Pathways
- Color-code genes based on stat results
- Integrated & interactive KEGG database for human, mouse, & rat as well as over 1,500 other organisms




Partek Pathway Extends Biological Interpretation

Dedicated, Intuitive Workflows for Ease of Use

Workflows **Gene Expression**


Gene Expression

- ▼ Import
 - Import samples ✓
 - Add sample attributes ✓
 - Edit sample information ✓
 - Choose sample ID column ✓
- ▶ QA/QC
- ▼ Analysis
 - Detect differentially expressed genes
 - Plot sources of variation
 - Create gene list
- ▼ Visualization
 - Cluster based on significant genes
 - Plot chromosome view
- ▼ Biological Interpretation
 - Gene set analysis
 - Pathway analysis
- ▶ Genomic Integration
- ▶ MicroRNA Integration



Exploratory Analysis

- High dimensional visualization
 - Principal components analysis (PCA)
 - Multidimensional scaling
- Clustering analysis
 - Hierarchical clustering
 - K-means clustering
 - Self organizing map
- Pattern searching
- Time series data analysis



Inferential Statistics

Parametric


- t-Test
- z-Test
- ANOVA
- Welch's ANOVA
- Pearson correlation

Non-Parametric


- Mann-Whitney
- Kruskal-Wallis
- Friedman
- Chi-square
- Rank correlation

Others

- Multiple test corrections
- Descriptive statistics



Survival Analysis



Kaplan-Meier curve

Categorical variable correlation with outcome (e.g. p53 status)

Cox regression

Continuous variable correlation with outcome (e.g. gene expression)

Calculate hazard ratio (HR)

Biomarker evaluation

Summary

	BaseSpace	Flow	Genomics Suite
Install	illumina 網頁版雲端分析	Partek 網頁版雲端分析	Partek 單機版分析
可分析項目	Only NGS data 超過80種不同應用的 分析流程(app)	Only NGS data DNA-seq, RNA-seq, ChIP-seq, Metagenomics	NGS, Microarray, qPCR Gene Expression, EPIC array, RNA- seq, ChIP-seq...etc.
特色	<ol style="list-style-type: none">1. 使用簡單2. 自動化分析流程3. 不需調整分析參數	<ol style="list-style-type: none">1. 使用簡單2. 較多統計及繪圖功能3. 可自建分析流程	<ol style="list-style-type: none">1.較多統計及繪圖功能2. 介面清楚易懂

THANK YOU

A decorative illustration of a branch with red and pink leaves and green berries, framing the text 'THANK YOU'. The text is written in a black, serif font. The word 'THANK' is on the top line and 'YOU' is on the bottom line. The branch starts on the left, goes up and right, then down and right, and finally down and left. It has several leaves in shades of red and pink, some with white vein patterns. There are also small green clusters of berries. The entire illustration is set against a plain white background.